


ターゲットバージョン指定/解除 設定手順 (JP)



手順 (1)


バッチファイルを右クリックし、表示されるコンテキストメニューからRun as Administratorを選択する。使用するバッチファイル名は以下の通りです。

ターゲットバージョン指定用：

SetTargeFeatureUpdateVersion_21H2.bat

ターゲットバージョン解除用：

RemoveTargeFeatureUpdateVersion.bat




手順 (2)

ダイアログ上で“SUCCESS”と表示されることを確認する。確認後、OKボタンでコマンドプロンプトを閉じて終了させる。

(注 1)

本バッチファイルと共通の設定をする [Local Group Policy Editor]内の[Select the target Feature Update version]設定がありますが、本バッチファイルで指定・解除を実行しても表示上の設定は変更されません。

Set / Cancel Target version setting procedure (EN)



Step 1

Right-click the batch file and select “Run as Administrator” from the displayed context menu. The batch file name to be used is as follows.

For target version setting:

SetTargeFeatureUpdateVersion_21H2.bat

For canceling the target version setting:

RemoveTargeFeatureUpdateVersion.bat

```
Administrator: C:\windows\System32\cmd.exe
2022/02/22 02:26:25, INFO, TargetReleaseVersion: 1
2022/02/22 02:26:25, INFO, TargetReleaseVersionInfo: 21H2
2022/02/22 02:26:25, INFO, DeferQualityUpdates: 1
2022/02/22 02:26:25, INFO, DeferQualityUpdatesPeriodInDays: 30
2022/02/22 02:26:25, INFO, PauseQualityUpdatesStartTime:
2022/02/22 02:26:25, DEBUG, GetRegSettings End
2022/02/22 02:26:25, INFO, -----Set the new configuration-----
2022/02/22 02:26:25, DEBUG, SetRegSettings Start
2022/02/22 02:26:25, DEBUG, GetRegSettings Start
2022/02/22 02:26:25, INFO, ProductVersion: Windows 10
2022/02/22 02:26:25, INFO, TargetReleaseVersion: 1
2022/02/22 02:26:25, INFO, TargetReleaseVersionInfo: 21H2
2022/02/22 02:26:25, INFO, DeferQualityUpdates: 1
2022/02/22 02:26:25, INFO, DeferQualityUpdatesPeriodInDays: 30
2022/02/22 02:26:25, INFO, PauseQualityUpdatesStartTime:
2022/02/22 02:26:25, DEBUG, GetRegSettings End
2022/02/22 02:26:25, DEBUG, SetRegSettings End
2022/02/22 02:26:25, INFO, -----Get the configuration After the change
2022/02/22 02:26:25, DEBUG, GetRegSettings Start
2022/02/22 02:26:25, INFO, ProductVersion: Windows 10
2022/02/22 02:26:25, INFO, TargetReleaseVersion: 1
2022/02/22 02:26:25, INFO, TargetReleaseVersionInfo: 21H2
2022/02/22 02:26:25, INFO, DeferQualityUpdates: 1
2022/02/22 02:26:25, INFO, DeferQualityUpdatesPeriodInDays: 30
2022/02/22 02:26:25, INFO, PauseQualityUpdatesStartTime:
2022/02/22 02:26:25, DEBUG, GetRegSettings End
2022/02/22 02:26:25, INFO, -----function is success-----
2022/02/22 02:26:25, INFO, Set TargetFeatureUpdateVersion is SUCCESS.
```

Step 2

Confirm that "SUCCESS" is displayed on the dialog. After confirmation, close the command prompt with the OK button and exit.

Note:

There is a [Select the target Feature Update version] setting in the [Local Group Policy Editor] that makes common settings with this batch file, but the settings on the display will not be changed even if you set / cancel in this batch file.