



## LMPLFLN5XBD

Product name	LMPLFLN5XBD
Product number	N2183600
UPC code	1-U2N342

Image



Dimension



Inquiry

[Contact us](#)

Product number	N2183600
----------------	----------

UPC code	1-U2N342
----------	----------

Major application	• Industrial
-------------------	--------------

Magnification [X]	5
Numerical aperture (NA)	0.13
Working distance (WD) [mm]	15
Objective field number [mm]	26.5
Cover glass thickness [mm]	-
Immersion medium	Air/Dry
Spring loaded	N/A
Correction collar	N/A

Remarks (Correction collar)

Iris	N/A
Correction level of chromatic aberration	Semiapochromat (FL)
Parfocalizing distance [mm]	45
Back focal plane (BFP) position [mm]	-3.0
Type of screw thread	W26×0.706
Brightfield (Reflected)	Good
Brightfield (Transmitted)	Good
Darkfield (Reflected)	Good
Darkfield (Transmitted)	N/A

Last updated: Jul 29 2024

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

## LMPLFLN5XBD

DIC (Reflected)	Good
DIC (Transmitted)	N/A
Phase contrast	N/A
Relief contrast	N/A
Polarization	Limitation
Fluorescence (B, G excitation)	N/A
UV fluorescence (at 365 nm)	N/A
Multiphoton	N/A
TIRF	N/A
IR	N/A
WLI	N/A
Auto focus	Yes