



LCACHN20XIPC

Product name	LCACHN20XIPC
Product number	N5229300
UPC code	1-U2PC145
Image	
Dimension	
Inquiry	Contact us
Web Store	Buy Now
Product number	N5229300
UPC code	1-U2PC145
Major application	• Life Science
Magnification [X]	20
Numerical aperture (NA)	0.4
Working distance (WD) [mm]	3.2
Objective field number [mm]	22
Cover glass thickness [mm]	1
Immersion medium	Air/Dry
Spring loaded	N/A
Correction collar	N/A
Remarks (Correction collar)	
Iris	N/A
Correction level of chromatic aberration	Achromat (ACH)
Parfocalizing distance [mm]	45
Back focal plane (BFP) position [mm]	-
Type of screw thread	W20.32×0.706 (RMS)
Brightfield (Reflected)	N/A
Brightfield (Transmitted)	Good
Darkfield (Reflected)	N/A

Last updated: May 15 2024

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

LCACHN20XIPC

Darkfield (Transmitted)	Limitation
DIC (Reflected)	N/A
DIC (Transmitted)	N/A
Phase contrast	Excellent
Relief contrast	N/A
Polarization	Limitation
Fluorescence (B, G excitation)	Limitation
UV fluorescence (at 365 nm)	N/A
Multiphoton	N/A
TIRF	N/A
IR	N/A
WLI	N/A
Auto focus	N/A