



## UPLFLN10X2

Product name	UPLFLN10X2
Product number	N5702900
UPC code	1-U2B5242
Image	
Dimension	
Inquiry	<a href="#">Contact us</a>
Web Store	<a href="#">Buy Now</a>
Product number	N5702900
UPC code	1-U2B5242
Major application	<ul style="list-style-type: none"> <li>Life Science</li> </ul>
Magnification [X]	10
Numerical aperture (NA)	0.3
Working distance (WD) [mm]	10
Objective field number [mm]	26.5
Cover glass thickness [mm]	-
Immersion medium	Air/Dry
Spring loaded	N/A
Correction collar	N/A
Iris	N/A
Correction level of chromatic aberration	Semiapochromat (FL)
Parfocalizing distance [mm]	45
Back focal plane (BFP) position [mm]	-19.1
Type of screw thread	W20.32×0.706 (RMS)
Brightfield (Reflected)	N/A
Brightfield (Transmitted)	Excellent
Darkfield (Reflected)	N/A
Darkfield (Transmitted)	Good

Last updated: Aug 01 2023

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

## UPLFLN10X2

DIC (Reflected)	N/A
DIC (Transmitted)	Excellent
Phase contrast	N/A
Relief contrast	N/A
Polarization	Good
Fluorescence (B, G excitation)	Excellent
UV fluorescence (at 365 nm)	Excellent
Multiphoton	N/A
TIRF	N/A
IR	N/A
WLI	N/A
Auto focus	Yes