



UPLFLN4X

Product name	UPLFLN4X
Product number	N5702800
UPC code	1-U2B5222
Image	
Dimension	
Inquiry	Contact us
Web Store	Buy Now
Product number	N5702800
UPC code	1-U2B5222
Major application	• Life Science
Magnification [X]	4
Numerical aperture (NA)	0.13
Working distance (WD) [mm]	17
Objective field number [mm]	26.5
Cover glass thickness [mm]	-
Immersion medium	Air/Dry
Spring loaded	N/A
Correction collar	N/A
Iris	N/A
Correction level of chromatic aberration	Semiapochromat (FL)
Parfocalizing distance [mm]	45
Back focal plane (BFP) position [mm]	-1.3
Type of screw thread	W20.32×0.706 (RMS)
Brightfield (Reflected)	N/A
Brightfield (Transmitted)	Excellent
Darkfield (Reflected)	N/A
Darkfield (Transmitted)	N/A

Last updated: Mar 05 2024

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

UPLFLN4X

DIC (Reflected)	N/A
DIC (Transmitted)	N/A
Phase contrast	N/A
Relief contrast	N/A
Polarization	Good
Fluorescence (B, G excitation)	Excellent
UV fluorescence (at 365 nm)	Excellent
Multiphoton	N/A
TIRF	N/A
IR	N/A
WLI	N/A
Auto focus	N/A