



UPLFLN60XOI

Product name	UPLFLN60XOI
Product number	N5703300
UPC code	1-U2B5322
Image	
Dimension	
Inquiry	Contact us
Web Store	Buy Now
Product number	N5703300
UPC code	1-U2B5322
Major application	<ul style="list-style-type: none"> Life Science
Magnification [X]	60
Numerical aperture (NA)	0.65 - 1.25
Working distance (WD) [mm]	0.12
Objective field number [mm]	26.5
Cover glass thickness [mm]	0.17
Immersion medium	Normal Oil
Spring loaded	Yes
Correction collar	N/A
Iris	Yes
Correction level of chromatic aberration	Semiapochromat (FL)
Parfocalizing distance [mm]	45
Back focal plane (BFP) position [mm]	-19.1
Type of screw thread	W20.32×0.706 (RMS)
Brightfield (Reflected)	N/A
Brightfield (Transmitted)	Excellent
Darkfield (Reflected)	N/A
Darkfield (Transmitted)	Good

Last updated: Aug 01 2023

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

UPLFLN60XOI

DIC (Reflected)	N/A
DIC (Transmitted)	Excellent
Phase contrast	N/A
Relief contrast	N/A
Polarization	Good
Fluorescence (B, G excitation)	Excellent
UV fluorescence (at 365 nm)	Excellent
Multiphoton	N/A
TIRF	N/A
IR	N/A
WLI	N/A
Auto focus	Yes